

# UV LAMPS

Lamps 6W

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-6.LC	365 / 254	145 x 48	2 x 6	UV Fluorescence
VL-6.L	365	145 x 48	1 x 6	UV Fluorescence
VL-6.C	254	145 x 48	1 x 6	UV Fluorescence
VL-206.BL	365	No Filter	2 x 6	Black Light
VL-206.G	254	No Filter	2 x 6	Germicidal

# UV LAMPS

Lamps 15W

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-215.L	365	295 x 76	2 x 15	UV Fluorescence
VL-215-M	312	295 x 76	2 x 15	UV Fluorescence
VL-215-C	254	295 x 76	2 x 15	UV Fluorescence
VL-115.L	365	295 x 66	1 x 15	UV Fluorescence
VL-115-M	312	295 x 66	1 x 15	UV Fluorescence
VL-115-C	254	295 x 66	1 x 15	UV Fluorescence
VL-315-BL	365	No Filter	3 x 15	Black Light
VL-315-G	254	No Filter	3 x 15	Germicidal
VL-215-BL	365	No Filter	2 x 15	Black Light
VL-215-G	254	No Filter	2 x 15	Germicidal
VL-115-BL	365	No Filter	1 x 15	Black Light
VL-115-G	254	No Filter	1 x 15	Germicidal

# UV LAMPS

Lamps 30W

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-330.BL	365	No Filter	3 x 30	Black Light
VL-330.G	254	No Filter	3 x 30	Germicidal
VL-230.BL	365	No Filter	2 x 30	Black Light
VL-230.G	254	No Filter	2 x 30	Germicidal
VL-130.BL	365	No Filter	1 x 30	Black Light
VL-130.G	254	No Filter	1 x 30	Germicidal



# UV LAMPS

*Irradiation Surface For A Height Of About 2m*

Power (W)	1 Tube	2 Tubes	3 Tubes
15	1.5 x 2m	1.5 x 3m	1.5 x 4m
20	2m x 2m	2m x 3m	2m x 4m
30	2.5m x 2.5m	2.5m x 3.5m	2.5m x 4.5m
40	3m x 3m	3m x 4m	3m x 5m