



UV TABLES

Transilluminators 8W

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
ECX-F26.MX V1	365	210 x 260	5 x 8	UV SuperBright
ECX-F26.M V1	312	210 x 260	5 x 8	UV
ECX-F26.C V1	254	210 x 260	5 x 8	UV
ECX-F20.L V1	365	200 x 200	5 x 8	UV
ECX-F20.M V1	312	200 x 200	5 x 8	UV
ECX-F20.C V1	254	200 x 200	5 x 8	UV
ECX-F20.BLUE V1	-	200 x 200	-	Blue Light
ECX-F20.W V1	-	No Filter	5 x 8	White Light
ECX-F15.M	312	150 x 150	4 x 8	UV
ECX-F15.C	254	150 x 150	4 x 8	UV



UV TABLES

Transilluminators 8W - Multiband

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
TCP-26.LMX V1	365 / 312	210 x 260	4 x 8 / 5 x 8	UV SuperBright
TCP-26.LM V1	365 / 312	210 x 260	4 x 8 / 5 x 8	UV
TCP-26.LC V1	365 / 254	210 x 260	4 x 8 / 5 x 8	UV
TCP-26.MC V1	312 / 254	210 x 260	4 x 8 / 5 x 8	UV
TCP-20.LM V1	365 / 312	200 x 200	4 x 8 / 5 x 8	UV
TCP-20.LC V1	365 / 254	200 x 200	4 x 8 / 5 x 8	UV
TCP-20.MC V1	312 / 254	200 x 200	4 x 8 / 5 x 8	UV